


## Introducing the Pannoramic MIDI II

Sysmex Deutschland GmbH is glad to invite you to an introduction of the Pannoramic MIDI II whole slide scanner for brightfield and fluorescence:

- High-quality digitalization of complete histopathological slides
- User-friendly software with both, full flexibility and efficient batch mode
- Multilayer (Z-Stack) scanning for highest content in BF and FL

[Please join our webinar on Microsoft Teams, Thursday Sept. 23<sup>rd</sup>, 3 p.m.](#)



If you want to have a personal view on the MIDI II please contact Dr. Antonio Virgilio Failla ([a.failla@uke.de](mailto:a.failla@uke.de)) to register for an on-site visit at the UKE room 01.068 between Oct. 20<sup>th</sup> and Nov. 19<sup>th</sup>.

Please have a first look at the Sysmex website:  
<https://www.sysmex.de/n/produkte/details/pannорamic-midi-ii-rx.html>

